



E-Chrom Tech Co., Ltd.

TEL: +886-2-2511-1887

FAX: +886-2-2537-1870

EMAIL: service@echromtech.com

ADD: 10F-12, NO.20, Min Chuan W. Road, Taipei
104, Taiwan

WEB: www.echromtech.com

B11- Shaking Bath



Descriptions & Features

- A. Microchip thermal control that enables to heat immediately and maintain stable temperature.
- B. Two 4-digit panels separately display set temperature and actual temperature, easy and direct to operate.
- C. Automatic power failure of max temperature guarantees safe and reliable operation.
- D. Error correction function of temperature controller allows more accurate temperature control.
- E. Unrestrainedly adjustable shaking times and duration.
- F. Steeple roof design avoids condensed water dropping into test tubes.

Specifications

Model	Shaking Bath
Temperature Range	Room temperature + 5°C ~ 100°C
Temperature Increment	±0.1°C
Temperature Accuracy	±0.2°C
Max Temperature Protection	Set temperature + 10°C
Heating Power Consumption (W)	600W ~ 1200W
Time Set	10min ~ 99hr 59min
Vibrating frequency	40 ~ 150 times/min
Vibrating displacement	18 mm (horizontal)
Vibrating time	10 ~ 90min
Capacity (L)	20L
Internal Dimensions (W × L × D)	490mm × 295mm × 150mm
External Dimensions (W × L × D)	500mm × 300mm × 300mm



E-Chrom Tech Co., Ltd.

TEL: +886-2-2511-1887

FAX: +886-2-2537-1870

EMAIL: service@echromtech.com

ADD: 10F-12, NO.20, Min Chuan W. Road, Taipei
104, Taiwan

WEB: www.echromtech.com

Optional accessories



(A)



(B)



(C)



(D)

Dimensions of stainless rack	360mm × 275mm × 150mm (W × L × D)
Adjustable springs of each rack	17 × 12
Test tube tray	Holds three test tube racks
Test tube rack	Holds $\Phi 13.5 \times 60$, $\Phi 20 \times 40$, $\Phi 25 \times 24$, $\Phi 30 \times 24$
Erlenmeyer flask tray	Holds 50ml × 24, 100ml × 15, 200ml × 6, 250ml × 6, 500ml × 6